

[Analysis Information]

Product Name:	TGX-221
Catalog No.:	BCC1244
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:					
Solvent A:					
Solvent B:					
Gradient:	<table border="1"><tr><th>Time</th><th>A</th></tr><tr><td></td><td>Stop</td></tr></table>	Time	A		Stop
Time	A				
	Stop				
Flow rate:	0.5ml/min				
Wavelength:					
Volume:	15µl				

[Result]

