

[Analysis Information]

Product Name:	Beta-Apopicropodophyllin
Catalog No.:	BCC1388
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Athena C18, 3 μ m, 2.1mm $\times$ 100 mm														
Solvent A:	0.1% H3PO4 in 100% Acetonitrile														
Solvent B:	0.1% H3PO4 in 100% Water														
Gradient:	<table border="1"> <thead> <tr> <th>Time</th> <th>A</th> </tr> </thead> <tbody> <tr> <td>0.0min</td> <td>10%</td> </tr> <tr> <td>3.5min</td> <td>95%</td> </tr> <tr> <td>8.0min</td> <td>95%</td> </tr> <tr> <td>9.0min</td> <td>10%</td> </tr> <tr> <td>10.0min</td> <td>10%</td> </tr> <tr> <td>10.0min</td> <td>Stop</td> </tr> </tbody> </table>	Time	A	0.0min	10%	3.5min	95%	8.0min	95%	9.0min	10%	10.0min	10%	10.0min	Stop
Time	A														
0.0min	10%														
3.5min	95%														
8.0min	95%														
9.0min	10%														
10.0min	10%														
10.0min	Stop														
Flow rate:	0.5ml/min														
Wavelength:	220nm														
Volume:	15 μ l														

[Result]

