

[Analysis Information]

Product Name: IRAK-1-4 Inhibitor I
 Catalog No.: BCC1659
 Operator: Jianhong Lee
 Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm $\times$ 100 mm

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

Time	A
0.0min	10%
10.0min	95%
15.0min	10%
18.0min	10%
18.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 254nm

Volume: 15 μ l

[Result]

