

[Analysis Information]

Product Name: MRT67307
Catalog No.: BCC1779
Operator: Jianhong Lee
Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm \times 100 mm

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

| Time | A |
|---------|------|
| 0.0min | 0% |
| 2.0min | 0% |
| 6.0min | 95% |
| 9.0min | 95% |
| 10.0min | 0% |
| 10.0min | Stop |

Flow rate: 0.5ml/min

Wavelength: 254nm

Volume: 15 μ l

[Result]

