

### [ Analysis Information ]

Product Name: SCH900776 S-isomer  
 Catalog No.: BCC1936  
 Operator: Jianhong Lee  
 Batch No.: b01

### [ Analysis Condition ]

Column: Athena C18, 3  $\mu$ m, 2.1mm  $\times$  100 mm

Solvent A: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Acetonitrile

Solvent B: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Water

| Time   | A    |
|--------|------|
| 0.0min | 0%   |
| 2.0min | 0%   |
| 7.0min | 95%  |
| 10.0   | 95%  |
| 10.0   | Stop |

Flow rate: 0.5ml/min

Wavelength: 230nm

Volume: 15 $\mu$ l

### [ Result ]

