

### [ Analysis Information ]

Product Name:	Tasquinimod
Catalog No.:	BCC1987
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column:							
Solvent A:							
Solvent B:							
Gradient:	<table border="1"><tr><th>Time</th><th>A</th></tr><tr><td></td><td></td></tr><tr><td></td><td>Stop</td></tr></table>	Time	A				Stop
Time	A						
	Stop						
Flow rate:	0.5ml/min						
Wavelength:							
Volume:	15µl						

### [ Result ]

