

[Analysis Information]

Product Name: PF-04929113 (SNX-5422)
 Catalog No.: BCC2130
 Operator: Jianhong Lee
 Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm $\times$ 100 mm

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

Gradient:	Time	A
	0.0min	10%
	3.5min	95%
	8.0min	95%
	9.0min	10%
	10.0min	10%
	10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 250nm

Volume: 15 μ l

[Result]

