

[Analysis Information]

Product Name:	VX-680 (MK-0457, Tozasertib)
Catalog No.:	BCC2167
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Athena C18, 3 μ m, 2.1mm \times 100 mm	
Solvent A:	0.1% H ₃ PO ₄ in 100% Acetonitrile	
Solvent B:	0.1% H ₃ PO ₄ in 100% Water	
Gradient:	Time	A
	0.0min	0%
	2.0min	0%
	7.0min	95%
	10.0min	95%
	10.0min	Stop
Flow rate:	0.5ml/min	
Wavelength:	254nm	
Volume:	15 μ l	

[Result]

