

### [ Analysis Information ]

|               |                                     |
|---------------|-------------------------------------|
| Product Name: | Baricitinib (LY3009104, INCB028050) |
| Catalog No.:  | BCC2195                             |
| Operator:     | Jianhong Lee                        |
| Batch No.:    | b01                                 |

### [ Analysis Condition ]

|             |                                                          |      |
|-------------|----------------------------------------------------------|------|
| Column:     | Athena C18, 3 $\mu$ m, 2.1mm $\times$ 100 mm             |      |
| Solvent A:  | 0.1% H <sub>3</sub> PO <sub>4</sub> in 100% Acetonitrile |      |
| Solvent B:  | 0.1% H <sub>3</sub> PO <sub>4</sub> in 100% Water        |      |
| Gradient:   | Time                                                     | A    |
|             | 0.0min                                                   | 0%   |
|             | 2.0min                                                   | 0%   |
|             | 7.0min                                                   | 95%  |
|             | 10.0min                                                  | 95%  |
|             | 10.0min                                                  | Stop |
| Flow rate:  | 0.5ml/min                                                |      |
| Wavelength: | 230nm                                                    |      |
| Volume:     | 15 $\mu$ l                                               |      |

### [ Result ]

