

### [ Analysis Information ]

Product Name: Pomalidomide (CC-4047)  
 Catalog No.: BCC2246  
 Operator: Jianhong Lee  
 Batch No.: b01

### [ Analysis Condition ]

Column: Athena C18, 3  $\mu$ m, 2.1mm  $\times$  100 mm

Solvent A: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Acetonitrile

Solvent B: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Water

Time	A
0.0min	10%
3.5min	95%
8.0min	95%
9.0min	10%
10.0min	10%
10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 259nm

Volume: 15 $\mu$ l

### [ Result ]

