

[Analysis Information]

Product Name:	Topiramate
Catalog No.:	BCC2314
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm \times 100 mm

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

Gradient:	Time	A
	0.0min	0%
	2.0min	0%
	6.0min	95%
	9.0min	95%
	10.0min	0%
	10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 205nm

Volume: 15 μ l

[Result]

