

[Analysis Information]

Product Name:	Tenovin-6
Catalog No.:	BCC3667
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column: Athena C18 (3 μ m, 2.1 mm \times 100 mm)

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

Time	A
0.0min	10%
2.0min	65%
3.5min	95%
4.5min	10%
5.5min	10%
10.0min	10%
10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 254nm

Volume: 15 μ l

[Result]

