

### [ Analysis Information ]

Product Name: GNF-5837  
Catalog No.: BCC3668  
Operator: Jianhong Lee  
Batch No.: b01

### [ Analysis Condition ]

Column: Athena C18, 3  $\mu$ m, 2.1mm  $\times$  100 mm  
Solvent A: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Acetonitrile Water  
Solvent B:  
Gradient:

Time	A
0.0min	10%
3.5min	95%
8.0min	95%
9.0min	10%
10.0min	10%
10.0min	Stop

Flow rate: 0.5ml/min  
Wavelength: 254nm  
Volume: 15 $\mu$ l

### [ Result ]

