

[Analysis Information]

| | |
|---------------|-----------------|
| Product Name: | Propranolol HCl |
| Catalog No.: | BCC4336 |
| Operator: | Jianhong Lee |
| Batch No.: | b01 |

[Analysis Condition]

| Column: | | | | | |
|-------------|---|------|---|--|------|
| Solvent A: | | | | | |
| Solvent B: | | | | | |
| Gradient: | <table border="1"><tr><th>Time</th><th>A</th></tr><tr><td></td><td>Stop</td></tr></table> | Time | A | | Stop |
| Time | A | | | | |
| | Stop | | | | |
| Flow rate: | 0.5ml/min | | | | |
| Wavelength: | | | | | |
| Volume: | 15 μ l | | | | |

[Result]

