

[Analysis Information]

Product Name:	Sotalol
Catalog No.:	BCC4356
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:					
Solvent A:					
Solvent B:					
Gradient:	<table border="1" style="width: 100%; text-align: center;"> <tr> <td style="width: 20%;">Time</td> <td>A</td> </tr> <tr> <td></td> <td>Stop</td> </tr> </table>	Time	A		Stop
Time	A				
	Stop				
Flow rate:	0.5ml/min				
Wavelength:					
Volume:	15µl				

[Result]

