

[Analysis Information]

Product Name:	FT-207 (NSC 148958)
Catalog No.:	BCC4455
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:					
Solvent A:					
Solvent B:					
Gradient:	<table border="1"><tr><td>Time</td><td>A</td></tr><tr><td></td><td>Stop</td></tr></table>	Time	A		Stop
Time	A				
	Stop				
Flow rate:	0.5ml/min				
Wavelength:					
Volume:	15 μ l				

[Result]

