

### [ Analysis Information ]

Product Name:	Donepezil HCl
Catalog No.:	BCC4569
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column:	Athena C18, 3 $\mu$ m, 2.1mm $\times$ 100 mm												
Solvent A:	0.1% H <sub>3</sub> PO <sub>4</sub> in 100% Acetonitrile												
Solvent B:	0.1% H <sub>3</sub> PO <sub>4</sub> in 100% Water												
Gradient:	<table border="1"><thead><tr><th>Time</th><th>A</th></tr></thead><tbody><tr><td>0.0min</td><td>0%</td></tr><tr><td>2.0min</td><td>0%</td></tr><tr><td>7.0min</td><td>95%</td></tr><tr><td>10.0min</td><td>95%</td></tr><tr><td>10.0min</td><td>Stop</td></tr></tbody></table>	Time	A	0.0min	0%	2.0min	0%	7.0min	95%	10.0min	95%	10.0min	Stop
Time	A												
0.0min	0%												
2.0min	0%												
7.0min	95%												
10.0min	95%												
10.0min	Stop												
Flow rate:	0.5ml/min												
Wavelength:	210nm												
Volume:	15 $\mu$ l												

### [ Result ]

