

[Analysis Information]

| | |
|---------------|---------------------------------|
| Product Name: | NMDA (N-Methyl-D-aspartic acid) |
| Catalog No.: | BCC4590 |
| Operator: | Jianhong Lee |
| Batch No.: | b01 |

[Analysis Condition]

| | | |
|-------------|--|------|
| Column: | Athena C18, 3 μ m, 2.1mm \times 100 mm | |
| Solvent A: | 0.1% H3PO4 in 100% Acetonitrile | |
| Solvent B: | 0.1% H3PO4 in 100% Water | |
| Gradient: | Time | A |
| | 0.0min | 0% |
| | 2.0min | 0% |
| | 6.0min | 95% |
| | 9.0min | 95% |
| | 10.0min | 0% |
| | 10.0min | Stop |
| Flow rate: | 0.5ml/min | |
| Wavelength: | 205nm | |
| Volume: | 15 μ l | |

[Result]

