

### [ Analysis Information ]

Product Name: Salicylanilide  
Catalog No.: BCC4712  
Operator: Jianhong Lee  
Batch No.: b01

### [ Analysis Condition ]

Column: Athena C18, 3  $\mu$ m, 2.1mm  $\times$  100 mm

Solvent A: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Acetonitrile

Solvent B: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Water

Time	A
0.0min	0%
2.0min	0%
7.0min	95%
10.0min	95%
10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 280nm

Volume: 15 $\mu$ l

### [ Result ]

