

[Analysis Information]

Product Name: Sodium Picosulfate
Catalog No.: BCC4845
Operator: Jianhong Lee
Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm \times 100 mm
Solvent A: 0.1% H3PO4 in 100% Acetonitrile
Solvent B: 0.1% H3PO4 in 100% Water

Gradient:	Time	A
	0.0min	0%
	2.0min	0%
	6.0min	95%
	9.0min	95%
	10.0min	0%
	10.0min	Stop

Flow rate: 0.5ml/min
Wavelength: 263nm
Volume: 15 μ l

[Result]

