

### [ Analysis Information ]

Product Name: Pirfenidone  
 Catalog No.: BCC5086  
 Operator: Jianhong Lee  
 Batch No.: b01

### [ Analysis Condition ]

Column: Athena C18, 3  $\mu$ m, 2.1mm  $\times$  100 mm

Solvent A: 0.1% H3PO4 in 100% Acetonitrile

Solvent B: 0.1% H3PO4 in 100% Water

Gradient:	Time	A
	0.0min	10%
	3.5min	95%
	8.0min	95%
	9.0min	10%
	10.0min	10%
	10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 315nm

Volume: 15 $\mu$ l

### [ Result ]

