

[Analysis Information]

Product Name: Beta-Lapachone
 Catalog No.: BCC5088
 Operator: Jianhong Lee
 Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm $\times$ 100 mm

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

Gradient:	Time	A
	0.0min	10%
	3.5min	95%
	8.0min	95%
	9.0min	10%
	10.0min	10%
	10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 254nm

Volume: 15 μ l

[Result]

