

[Analysis Information]

Product Name: Talabostat mesylate
 Catalog No.: BCC5357
 Operator: Jianhong Lee
 Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm $\times$ 100 mm

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

Time	A
0.0min	0%
2.0min	0%
6.0min	95%
9.0min	95%
10.0min	0%
10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 205nm

Volume: 15 μ l

[Result]

