

### [ Analysis Information ]

Product Name:	Meclofenamate Sodium
Catalog No.:	BCC5490
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column:	Athena C18, 3 $\mu$ m, 2.1mm $\times$ 100 mm														
Solvent A:	0.1% H <sub>3</sub> PO <sub>4</sub> in 100% Acetonitrile														
Solvent B:	0.1% H <sub>3</sub> PO <sub>4</sub> in 100% Water														
Gradient:	<table border="1"><thead><tr><th>Time</th><th>A</th></tr></thead><tbody><tr><td>0.0min</td><td>10%</td></tr><tr><td>3.5min</td><td>95%</td></tr><tr><td>8.0min</td><td>95%</td></tr><tr><td>9.0min</td><td>10%</td></tr><tr><td>10.0min</td><td>10%</td></tr><tr><td>10.0min</td><td>Stop</td></tr></tbody></table>	Time	A	0.0min	10%	3.5min	95%	8.0min	95%	9.0min	10%	10.0min	10%	10.0min	Stop
Time	A														
0.0min	10%														
3.5min	95%														
8.0min	95%														
9.0min	10%														
10.0min	10%														
10.0min	Stop														
Flow rate:	0.5ml/min														
Wavelength:	240nm														
Volume:	15 $\mu$ l														

### [ Result ]

