

### [ Analysis Information ]

Product Name: PFI-2  
Catalog No.: BCC5561  
Operator: Jianhong Lee  
Batch No.: b01

### [ Analysis Condition ]

Column: Athena C18, 3  $\mu$ m, 2.1mm  $\times$  100 mm  
Solvent A: 0.1% H3PO4 in 100% Acetonitrile  
Solvent B: 0.1% H3PO4 in 100% Water

Gradient:	Time	A
	0.0min	0%
	2.0min	0%
	6.0min	95%
	9.0min	95%
	10.0min	0%
	10.0min	Stop

Flow rate: 0.5ml/min  
Wavelength: 230nm  
Volume: 15 $\mu$ l

### [ Result ]

