

[Analysis Information]

Product Name: J 147
 Catalog No.: BCC6360
 Operator: Jianhong Lee
 Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm $\times$ 100 mm
 Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile
 Solvent B: 0.1% H₃PO₄ in 100% Water
 Gradient:

Time	A
0.0min	5%
9.0min	95%
15.0min	95%
15.0min	Stop

 Flow rate: 0.5ml/min
 Wavelength: 280nm
 Volume: 15 μ l

[Result]

