

### [ Analysis Information ]

Product Name: AZ5104  
 Catalog No.: BCC6389  
 Operator: Jianhong Lee  
 Batch No.: b01

### [ Analysis Condition ]

Column: Athena C18, 3  $\mu$ m, 2.1mm  $\times$  100 mm

Solvent A: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Acetonitrile

Solvent B: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Water

| Gradient: | Time    | A    |
|-----------|---------|------|
|           | 0.0min  | 0%   |
|           | 2.0min  | 0%   |
|           | 6.0min  | 95%  |
|           | 9.0min  | 95%  |
|           | 10.0min | 0%   |
|           | 10.0min | Stop |

Flow rate: 0.5ml/min

Wavelength: 271nm

Volume: 15 $\mu$ l

### [ Result ]

