


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	Puromycin dihydrochloride
<b>Catalog Number:</b>	BCC7860
<b>Cas No.:</b>	58-58-2
<b>Formula:</b>	C <sub>22</sub> H <sub>31</sub> Cl <sub>2</sub> N <sub>7</sub> O <sub>5</sub>
<b>M.Wt:</b>	544.43
<b>Purity:</b>	>98%
<b>Physical Description:</b>	N/A
<b>Synonyms:</b>	CL13900 dihydrochloride
<b>Chemical Name:</b>	3'-[α-Amino- <i>p</i> -methoxyhydrocinnamamido]-3'-deoxy- <i>N,N</i> -dimethyladenosine dihydrochloride
<b>SMILES:</b>	<chem>CN(C)C1=NC=NC2=C1N=CN2C3C(C(C(O3)CO)NC(=O)C(CC4=CC=C(C(=C4)OC)N)O.Cl.Cl</chem>
<b>Solubility:</b>	DMSO : ≥ 77.5 mg/mL (142.35 mM) H <sub>2</sub> O : 50 mg/mL (91.84 mM; Need ultrasonic and warming) *"≥" means soluble, but saturation unknown.



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure*

*that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893