


---

## PHYSICAL AND CHEMICAL PROPERTIES

---

<b>Product Name:</b>	6-Angeloyloxyditropan-3-yl itaconate
<b>Catalog Number:</b>	BCN1867
<b>Cas No.:</b>	182015-05-0
<b>Formula:</b>	C <sub>26</sub> H <sub>38</sub> N <sub>2</sub> O <sub>6</sub>
<b>M.Wt:</b>	474.60
<b>Purity:</b>	>98%
<b>Physical Description:</b>	Oil
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	1-O-[(1S,5R)-8-methyl-8-azabicyclo[3.2.1]octan-3-yl] 4-O-[(1R,3R,5S,6R)-8-methyl-6-[(Z)-2-methylbut-2-enoyl]oxy-8-azabicyclo[3.2.1]octan-3-yl] 2-methylidenebutanedioate
<b>SMILES:</b>	<chem>CC=C(C)C(=O)OC1CC2CC(CC1N2C)OC(=O)CC(=C)C(=O)OC3CC4CCC(C3)N4C</chem>
<b>Solubility:</b>	N/A



biocrick.com

---

## STORAGE AND SHIPPING

---

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

---

## INTENDED USE

---

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

---

## CAUTION!

---

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893