


## PHYSICAL AND CHEMICAL PROPERTIES

|                              |  |
|------------------------------|--|
| <b>Product Name:</b>         | 6-Methoxykaempferol 3-O-rutinoside   |
| <b>Catalog Number:</b>       | BCN2379  |
| <b>Cas No.:</b>              | 403861-33-6  |
| <b>Formula:</b>              | C <sub>28</sub> H <sub>32</sub> O <sub>16</sub>  |
| <b>M.Wt:</b>                 | 624.6  |
| <b>Purity:</b>               | >98%   |
| <b>Physical Description:</b> | Yellow powder  |
| <b>Synonyms:</b>             | N/A  |
| <b>Chemical Name:</b>        | 5,7-dihydroxy-2-(4-hydroxyphenyl)-6-methoxy-3-[(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-[(2R,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl]oxymethyl]oxan-2-yl]oxychromen-4-one |
| <b>SMILES:</b>               | <chem>CC1C(C(C(C(O1)OCC2C(C(C(C(O2)OC3=C(OC4=CC(=C(C(=C4C3=O)O)OC)O)C5=CC=C(C(=C5)O)O)O)O)O)O)O</chem>   |
| <b>Solubility:</b>           | N/A  |



## STORAGE AND SHIPPING

|                            |   |
|----------------------------|---|
| <b>Storage:</b>            | Store the product in sealed, cool and dry condition.  |
| <b>General tips:</b>       | For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months. |
| <b>Shipping Condition:</b> | Ship with RT , or blue ice upon request   |

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15,Building B2, Tianfu Life Science Park,88 South Keyuan Road, Hi-Tech Zone,Chengdu,Sichuan 610041,PRC  
Tech Support:service@biocrick.com Tel:+86-28-8543-3893 FAX:+86-28-8543-3893