

[Analysis Information]

Product Name: (-)-Usnic acid
Catalog No.: BCN4306
Operator: Jianhong Lee
Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm \times 100 mm
Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile
Solvent B: 0.1% H₃PO₄ in 100% Water
Gradient:

Time	A
0.0min	0%
2.0min	0%
6.0min	95%
9.0min	95%
10.0min	0%
10.0min	Stop

Flow rate: 0.5ml/min
Wavelength: 220nm
Volume: 15 μ l

[Result]

