


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	Hyperoside
<b>Catalog Number:</b>	BCN5570
<b>Cas No.:</b>	482-36-0
<b>Formula:</b>	C <sub>21</sub> H <sub>20</sub> O <sub>12</sub>
<b>M.Wt:</b>	464.4
<b>Purity:</b>	>98%
<b>Physical Description:</b>	Light yellow powder
<b>Synonyms:</b>	Hyperin; 3,3',4',5,7-Pentahydroxyflavone 3-galactoside; Quercetin 3-galactoside
<b>Chemical Name:</b>	2-(3,4-dihydroxyphenyl)-5,7-dihydroxy-3-[(2S,3R,4S,5R,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxychromen-4-one
<b>SMILES:</b>	<chem>O=C1C(=C(C=C1C2=C(C(=O)C3=C(C=C(C=C3O2)O)O)OC4C(C(C(C(O4)CO)O)O)O)O)O</chem>
<b>Solubility:</b>	DMSO : 250 mg/mL (538.35 mM; Need ultrasonic)



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15,Building B2, Tianfu Life Science Park,88 South Keyuan Road, Hi-Tech Zone,Chengdu,Sichuan 610041,PRC  
Tech Support:service@biocrick.com Tel:+86-28-8543-3893 FAX:+86-28-8543-3893