


## PHYSICAL AND CHEMICAL PROPERTIES

|                              |   |
|------------------------------|---|
| <b>Product Name:</b>         | Pachymic acid   |
| <b>Catalog Number:</b>       | BCN6347   |
| <b>Cas No.:</b>              | 29070-92-6  |
| <b>Formula:</b>              | C <sub>33</sub> H <sub>52</sub> O <sub>5</sub>  |
| <b>M.Wt:</b>                 | 528.76  |
| <b>Purity:</b>               | >98%  |
| <b>Physical Description:</b> | White powder  |
| <b>Synonyms:</b>             | 3-O-Acetyltumulosic acid  |
| <b>Chemical Name:</b>        | 2-(3-acetyloxy-16-hydroxy-4,4,10,13,14-pentamethyl-2,3,5,6,7,11,12,15,16,17-decahydro-1H-cyclopenta[a]phenanthren-17-yl)-6-methyl-5-methylideneheptanoic acid |
| <b>SMILES:</b>               | <chem>CC(C)C(=C)CCC(C1C(CC2(C1(CCC3=C2CCC4C3(CCC(C4(C)C)OC(=O)C)C)C)C)O)C(=O)O</chem>   |
| <b>Solubility:</b>           | DMSO : 7.14 mg/mL (13.50 mM; Need ultrasonic)<br>H <sub>2</sub> O : < 0.1 mg/mL (insoluble)   |



## STORAGE AND SHIPPING

|                            |   |
|----------------------------|---|
| <b>Storage:</b>            | Store the product in sealed, cool and dry condition.  |
| <b>General tips:</b>       | For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months. |
| <b>Shipping Condition:</b> | Ship with RT , or blue ice upon request   |

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure*

*that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893