


---

## PHYSICAL AND CHEMICAL PROPERTIES

---

<b>Product Name:</b>	Glucogenin C mono-D-thevetoside
<b>Catalog Number:</b>	BCN7089
<b>Cas No.:</b>	849201-84-9
<b>Formula:</b>	C <sub>28</sub> H <sub>40</sub> O <sub>9</sub>
<b>M.Wt:</b>	520.61
<b>Purity:</b>	>98%
<b>Physical Description:</b>	Powder
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	(4S,5R,8S,13R,16S,19R,22R)-8-[(2R,3R,4S,5R,6R)-3,5-dihydroxy-4-methoxy-6-methyloxan-2-yl]oxy-5,19-dimethyl-15,18,20-trioxapentacyclo[14.5.1.04,13.05,10.019,22]docosa-1(21),10-dien-14-one
<b>SMILES:</b>	<chem>CC1C(C(C(C(O1)OC2CCC3(C4CCC5=COC6(C5C(CO6)OC(=O)C4CC=C3C2)C)C)O)OC)O</chem>
<b>Solubility:</b>	N/A



---

## STORAGE AND SHIPPING

---

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

---

## INTENDED USE

---

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

---

## CAUTION!

---

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893