


---

## PHYSICAL AND CHEMICAL PROPERTIES

---

|                              |   |
|------------------------------|---|
| <b>Product Name:</b>         | Ganoderenic acid E  |
| <b>Catalog Number:</b>       | BCN8241   |
| <b>Cas No.:</b>              | 110241-23-1   |
| <b>Formula:</b>              | C <sub>30</sub> H <sub>40</sub> O <sub>8</sub>  |
| <b>M.Wt:</b>                 | 528.63  |
| <b>Purity:</b>               | >98%  |
| <b>Physical Description:</b> | Powder  |
| <b>Synonyms:</b>             | N/A   |
| <b>Chemical Name:</b>        | (E,2S)-6-[(7S,10S,12S,13R,14R)-7,12-dihydroxy-4,4,10,13,14-pentamethyl-3,11,15-trioxo-1,2,5,6,7,12,16,17-octahydrocyclopenta[a]phenanthren-17-yl]-2-methyl-4-oxohept-5-enoic acid |
| <b>SMILES:</b>               | <chem>CC(CC(=O)C=C(C)C1CC(=O)C2(C1(C(C(=O)C3=C2C(CC4C3(CCC(=O)C4(C)C)C)O)O)C)C(=O)O</chem>  |
| <b>Solubility:</b>           | N/A   |



---

## STORAGE AND SHIPPING

---

|                            |   |
|----------------------------|---|
| <b>Storage:</b>            | Store the product in sealed, cool and dry condition.  |
| <b>General tips:</b>       | For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months. |
| <b>Shipping Condition:</b> | Ship with RT , or blue ice upon request   |

---

## INTENDED USE

---

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

---

## CAUTION!

---

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893