

## PHYSICAL AND CHEMICAL PROPERTIES

**Product Name:** Polyporusterone B

**Catalog Number:** BCX0680

**Cas No.:** 141360-89-6

**Formula:** C<sub>28</sub>H<sub>44</sub>O<sub>6</sub>

**M.Wt:** 476.65

**Purity:** >98%


**Physical Description:** N/A

**Synonyms:** N/A

**Chemical Name:** (2S,3R,5R,9R,10R,13R,14S,17S)-17-[(2R,3R)-2,3-dihydroxy-6-methyl-5-methylideneheptan-2-yl]-2,3,14-trihydroxy-10,13-dimethyl-2,3,4,5,9,11,12,15,16,17-decahydro-1H-cyclopenta[a]phenanthren-6-one

**SMILES:** CC(C)C(=C)CC(C(C)(C1CCC2(C1(CCC3C2=CC(=O)C4C3(CC(C(C4)O)O)C)O)O)O)O

**Solubility:** N/A



## STORAGE AND SHIPPING

**Storage:** Store the product in sealed, cool and dry condition.

**General tips:** For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.

**Shipping Condition:** Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15,Building B2, Tianfu Life Science Park,88 South Keyuan Road, Hi-Tech Zone,Chengdu,Sichuan 610041,PRC  
Tech Support:service@biocrick.com Tel:+86-28-8543-3893 FAX:+86-28-8543-3893